

Probe generation by in vitro transcription (RIBOPROBES)

Probes can be generated by in vitro transcription reaction in the presence of labeled nucleotides. The target DNA fragment to be detected is cloned into the polylinker region of a transcription vector containing promoters for either SP6, T3 or T7 RNA polymerase. To synthesize an anti-sense “run-off” transcript, a specific restriction endonuclease should be used to linearize the DNA template and to create a 5’ overhang on the anti-sense strand before transcription. The procedure described here incorporates one biotinylated nucleotide for approximately for every 20-25 nucleotides in the newly-synthesized transcript.

Equipment and reagents

- 10% SDS
- Proteinase K (10 mg/ml) (Roche Molecular Biochemicals)
- Phenol/chloroform/isoamyl alcohol (25:24:1) (Gibco)
- SP6, or T3, or T7 polymerase, 20 units/ μ l (Roche Molecular Biochemicals)
- Ethanol
- Biotin-16-UTP 10 mM (Roche Molecular Biochemicals)
- NTPs 10 mM (Roche Molecular Biochemicals)

Method

1. Perform restriction digest of 20-25 μ g of template DNA with specific restriction endonuclease which cut the template ~200-300 nt downstream of the promoter region to create an 5’ overhang to linearize the DNA template used for “run-off” transcription in 1.5 ml test tube
2. After the restriction digest add 10 μ l 10% SDS and 2 μ l proteinase K (10 mg/ml) and add enough water to make a total reaction volume of 200 μ l
3. Incubate the tube for 30 min at 50°C
4. Add 8.5 μ l 0.5 M EDTA (pH 8.0), 1.5 μ l water and 210 μ l phenol/chloroform/isoamyl alcohol (25:24:1). Vortex for 1 min
5. Centrifuge the tube at 13,000 x g for 5 min at room temperature

6. Collect upper aqueous layer into new 1.5 ml test tube. Add 200 μ l phenol/chloroform/isoamyl alcohol (25:24:1). Vortex for 45 seconds
7. Centrifuge the tube at 13,000 x g for 5 min at room temperature
8. Collect 180 μ l upper aqueous layer. Add 20 μ l 3M sodium acetate (pH 5.2) and 500 μ l pre-chilled (-20°C) 100% ethanol. Let the precipitate form for at least 20 min at -20°C
9. Centrifuge the tube at 13,000 x g for 30 min at 4°C
10. Discard the supernatant. Add 200 μ l pre-chilled (-20°C) 70% ethanol
11. Centrifuge the tube at 13,000 x g for 30 min at 4°C
12. Discard the supernatant
13. Dry the pellet under vacuum. Dissolve the pellet into 10 μ l nuclease-free water. Adjust the final concentration to $\sim 1 \mu\text{g}/\mu\text{l}$. If desired, store at -20°C
14. To set up in vitro transcription reaction add the following to a 1.5 ml microcentrifuge tube on ice, 2 μ l 10x transcription buffer (supplied with SP6, T3 or T7 RNA polymerase), 2 μ l 10 mM ATP, 2 μ l 10 mM CTP, 2 μ l 10 mM GTP, 1.3 μ l 10 mM UTP, 0.7 μ l 10 mM biotin-16-UTP, 1 μ g purified linearized DNA template, 0.5 μ l RNA pol (SP6, T3 or T7) (20U/ μ l). Enough nuclease-free water to make a total reaction volume of 20 μ l.
15. Mix the components and centrifuge briefly.
16. Incubate the tube for 2 h at 37°C
17. To the reaction tube, add 2 μ l 0.2 M EDTA (pH 8.0) to stop the reaction, 2 μ l 4M LiCl and 60 μ l pre-chilled (-20°C) 100% ethanol. Let the precipitate form for at least 20 min at -20°C
18. Centrifuge the tube at 13,000 x g for 15 min at 4°C
19. Discard the supernatant
20. Wash the pellet with 50 μ l pre-chilled 70% (v/v) ethanol
21. Centrifuge the tube at 13,000 x g for 10 min
22. Discard the supernatant
23. Dry the pellet under vacuum. Dissolve the RNA pellet for 30 min at 37°C in 50 μ l nuclease-free water plus 1 μ l RNase inhibitor (20U/ μ l)

24. To estimate the yield of the transcript. Run an aliquot of the transcript on an agarose gel beside an RNA standard of known concentration. Stain with ethidium bromide. Compare the relative intensity of staining between the labeled transcripts and the known standard^b.

^aIf you want to remove the template DNA, add 2 units DNase I, RNase-free to the test tube and incubate for 15 min at 37°C and then go to the step 4

^bIf you are not going to use the labeled probe immediately, store the probe solution at -70°C