

***Preparation of DNA probes for chromosome FISH: Whole chromosome painting probes labeling by DOP-PCR***

***Reagents and equipment***

- Source DNA (flow sorted or microdisected chromosomes)
- PCR buffer 10X Perkin Elmer without  $MgCl_2$  (Roche Molecular Biochemical)
- $MgCl_2$  solution 25 mM Perkin Elmer (Roche Molecular Biochemical)
- Taq DNA polymerase 5U/ $\mu$ l Perkin Elmer (Roche Molecular Biochemical)
- Deoxy-nucleotides dATP, dCTP, dGTP, and dTTP 100 mM (Roche Molecular Biochemical)
- Stock dNTPs solution for labeling PCR contains 0.2 mM of dATP, dCTP, dGTP and 0.15 mM of dTTP
- 0.06mM Fluorescein–dUTP. (Boehringer Mannheim) or Biotin-16-dUTP (Boehringer Mannheim) or Spectrum red-dUTP (Vysis)
- Universal primer for human genomic DNA amplification: UN1 (Midland Certified Reagent Co) Telenius [5'-CCGACTCGAGNNNNNNATGTGG-3'] or universal primer for mouse genomic DNA amplification: 22-mer (Midland Certified Reagent Co.) [5'-CGG ACT CGA GNN NNN NTA CAC C-3']
- Hi-Lo DNA marker (Minnesota Molecular)
- Tris acetate buffer
- Sodium Acetate 3M pH 5.2
- PCR thermocycler
- Microcentrifuge
- Microcentrifuge tubes

***Method***

1. Set up the PCR reaction.

PCR reaction mix:

| Reagent | Quantity ( $\mu$ l) |
|---------|---------------------|
|---------|---------------------|

|                                   |    |
|-----------------------------------|----|
| PCR buffer 10X                    | 10 |
| MgCl <sub>2</sub> 25mM            | 8  |
| Stock dNTPs solution <sup>a</sup> | 5  |
| Biotin–dUTP <sup>a</sup>          | 5  |
| dH <sub>2</sub> O                 | 65 |
| DNA (100-150 ng/μl)               | 4  |
| Primer (100 μM)                   | 2  |
| Taq polymerase (5U/μl)            | 1  |

## 2. Run the PCR reaction

| Step | Temperature (°C)     | Time (min)                |
|------|----------------------|---------------------------|
| 1    | 94                   | 1                         |
| 2    | 56                   | 1                         |
| 3    | 72                   | 3                         |
|      |                      | (+1 additional sec/cycle) |
| 4    | Steps 1-3 (29 times) |                           |
| 5    | 72                   | 10                        |

After the final step hold the PCR samples at 4°C until they are used.

3. To analyze the DOP-PCR products, mix 8 μl of the reaction products with 2 μl agarose gel loading buffer.
4. Apply the sample to a 1% w/v agarose gel in 1x TAE buffer.
5. Apply 10 μl of Hi-Lo DNA marker.
6. Run the gel for 45 min at 70 V/cm in 1x TAE buffer.
7. Stain the gel with ethidium bromide and observe in UV transilluminator. In the reaction samples you should see a smear of DNA ranging from about 200–500 bp.
8. Add to the remnant of the DOP-PCR labeled DNA 1/10 volume of 3M sodium acetate pH 5.2 and 3 volumes of cold absolute ethanol.
9. Put tubes at -70 °C for 30 min.

10. Spin down the samples at 14,000 g in 4 °C microcentrifuge for 30 min.
11. Remove carefully the supernatant and dry the DNA under vacuum for 3 min.
12. Resuspend the DNA in sterile water at a final concentration of 100 ng/μl.
13. Store the labeled DNA at –20 °C.

<sup>a</sup>Note that the concentration of the dNTPs and labeled-dUTP varies.

- Remnants of cytoplasm can impair the access of the probe to the target DNA of metaphase and especially interphase chromosomes. To improve the hybridization quality a treatment with 70% acetic acid and/or mild pepsin treatment is frequently performed. Here is the protocol

#### *Equipment and reagents*

- Glacial acetic acid
- Hydrochloric acid 1N
- Pepsin (Sigma)
- 1XPBS
- Coplin jars
- Slide warmer or waterbath

#### *Method for acetic acid treatment*

1. Soak the slide in an acetic acid 70 % v/v solution for 40-60 sec.
2. Rinse the slide in 1xPBS for 5 min at room temperature, shaking gently.
3. Dehydrate the slide through an ethanol series of 70 % v/v, 85% v/v, and 100 % v/v for 2 min each at room temperature and then air dry.

#### *Method for pepsin treatment*

1. Rinse the slide in 1XPBS for 5 min.

2. Apply to a 60 mm<sup>2</sup> coverslip 100 µl of a 10µg/ml solution of pepsin in 10 mM HCl. Touch the slide to the coverslip. Incubate in moist chamber at 37°C, for 5 min.
3. Rinse the slide in 1XPBS for 5 min.
4. Dehydrate the slide through an ethanol series of 70 % v/v, 85% v/v, and 100 % v/v for 2 min each at room temperature and then air dry.